

1. Alignment files

a. Upload

- i. Enable upload annotation files
- ii. FASTA file formats: .fna, .fa, .fasta (note: **.gz is not supported**)
- iii. Gene annotation file formats: .gtf, .gff, .gff3, .gtf.gz, .gff.gz, .gff3.gz
- iv. Wait for each file to finish uploading completely before uploading the next or proceeding.
(note: Especially for FASTA files which can be very big)




Fig 4.1: Upload incomplete




Fig 4.2: Upload complete

- v. Select correct file format for gene annotation file uploaded
(note: .gtf.gz → select “gtf” / .gff.gz → select “gff” / .gff3.gz → select “gff3”)
- vi. Click ‘Prepare Annotation Files’ button once done




Fig 4.3: Prepare annotation files

b. Download from NCBI

- i. In the NCBI tab, search the organism name and find the ncbi accession number.

NCBI Database Counts

Show 10 entries

Search for organism Search:

refseq_category	taxid	species_taxid	organism_name	infraspecific_name	isolate
All		All	Drosop	All	A

reference genome	7227	7227	Drosophila melanogaster		
------------------	------	------	-------------------------	--	--

Fig 4.4: Search organism

NCBI Database Counts

Show 10 entries

Search:

assembly_accession	bioproject	biosample	wgs_master	refseq_category	taxid	spec
All	A	All	All	All		All

1	GCF_000001215.4	PRJNA164	SAMN02803731	reference genome	7227	
---	-----------------	----------	--------------	------------------	------	--

Find accession number

Fig 4.5: Search accession no.

- ii. Enter ncbi accession number in the left panel

Upload Settings

Upload Annotation Files? ☐ Disable

Enter NCBI Accession No.

GCF_000001215.4

Fig 4.5: Enter accession no.

- iii. FASTA file: From the dropdown, select the option "XXX_genomic.fna.gz"

Select NCBI Fasta File

GCF_000001215.4_Release_6_plus_ISO1_MT_genomic.fna.gz

GCF_000001215.4_Release_6_plus_ISO1_MT_cds_from_genomic.fna.gz

GCF_000001215.4_Release_6_plus_ISO1_MT_genomic.fna.gz

GCF_000001215.4_Release_6_plus_ISO1_MT_protein.faa.gz

GCF_000001215.4_Release_6_plus_ISO1_MT_rna.fna.gz

GCF_000001215.4_Release_6_plus_ISO1_MT_rna_from_genomic.fna.gz

GCF_000001215.4_Release_6_plus_ISO1_MT_translated_cds.faa.gz

Fig 4.6: NCBI fasta file

- iv. Gene annotation file: From the dropdown, select the options “XXX_genomic.gtf.gz” (gene annotation format: “gtf”) or “XXX_genomic.gff.gz” (gene annotation format: “gff”)

Select NCBI Annotation File

GCF_000001215.4_Release_6_plus_ISO1_MT_genomic.gff.gz

Select Annotation File Format

☐ gtf ☒ gff

Fig 4.7: Find gff file

Select NCBI Annotation File

GCF_000001215.4_Release_6_plus_ISO1_MT_genomic.gtf.gz

Select Annotation File Format

☒ gtf ☐ gff

Fig 4.8: Find gtf file

- v. Click ‘Prepare Annotation Files’ button once done

Upload Settings

Upload Annotation Files? ☐

Enter NCBI Accession No.

GCF_000001215.4

Select NCBI Faasta File

GCF_000001215.4_Release_6_plus_ISO1_MT_genomic.fna.gz

Select NCBI Annotation File

GCF_000001215.4_Release_6_plus_ISO1_MT_genomic.gff.gz

Select Annotation File Format

☐ gtf ☒ gff ☐ gff3

[Prepare Annotation Files](#)

Fig 4.9: Prepare annotation files

- c. Sources of FASTA & gene annotation files
 - i. Ensembl (however files from ensembl take a long time to download so Ensembl download functionality is not supported in this shiny application)
 - ii. NCBI
- d. Importance of ‘Prepare Annotation’ functionality
 - i. Some chromosomes in fasta file not found in gtf file (note: Often when downloading _genomic.fna.gz from NCBI / dna.primary_assembly.fa.gz or dna.toplevel.fa.gz from Ensembl)

```
> length(names(readDNAStringSet('/Users/paigepattimusa/Desktop/shiny_1/r_4_1/
Drosophila_melanogaster.BDGP6.32.dna.toplevel.fa')))
[1] 1870
```

Fig 4.10: Fasta file from Ensembl

```
> length(seqlevels(makeTxDbFromGFF('/Users/paigepattimusa/Desktop/shiny_1/
r_4_1/Drosophila_melanogaster.BDGP6.32.107.gtf.gz')))
Import genomic features from the file as a GRanges object ... OK
Prepare the 'metadata' data frame ... OK
Make the TxDb object ... OK
[1] 29
```

Fig 4.11: GTF file from Ensembl

- ii. Some chromosomes in gtf file not found in FASTA file
(note: When downloading individual chromosome fasta files from Ensembl)

Index of /pub/release-107/fasta/mus_musculus/dna/

Individual fasta files		
04-Jun-2022 10:27	4339	
04-Jun-2022 08:49	58349084	
04-Jun-2022 08:51	38648477	
04-Jun-2022 08:49	36137719	
04-Jun-2022 08:49	35553040	
04-Jun-2022 08:49	35687980	
04-Jun-2022 08:49	36857618	
04-Jun-2022 08:51	30695959	
04-Jun-2022 08:51	28847171	
04-Jun-2022 08:49	27906831	
04-Jun-2022 08:49	26641990	
04-Jun-2022 08:49	17732438	
04-Jun-2022 08:49	54235200	
04-Jun-2022 08:49	47576831	
04-Jun-2022 08:49	46381857	
04-Jun-2022 08:49	44992076	
04-Jun-2022 08:51	44451585	
04-Jun-2022 08:49	42934232	
04-Jun-2022 08:49	38201500	
04-Jun-2022 08:51	36852135	
04-Jun-2022 08:49	5300	
04-Jun-2022 08:49	49575889	
04-Jun-2022 08:51	26759373	
04-Jun-2022 08:49	1394605	
04-Jun-2022 08:50	806418890	
04-Jun-2022 08:50	806418890	

Fig 4.12: Individual fasta files from Ensembl

- iii. Ensures that in these 2 instances, chromosomes in fasta file & gtf file are identical to prevent any errors later on when aligning.

2. In the “Settings” tab of the left side panel, select the sample for aligning.

Upload

Settings

Select FASTQ Sample

male_a

male_a

female_a

male_b

female_b

Fig 4.13: Select sample for aligning

3. Configure trim settings

(note: If spliced alignment is enabled and Rbowtie is chosen as the aligner, SpliceMap will be used which takes much longer than Rhisat2!)




Fig 4.14: Select aligner for spliced alignment




Fig 4.15: Use trimmed / untrimmed fastq files

Upload Settings

Select FASTQ Sample

male_a

Use Trimmed/Untrimmed Files?

☒ Trimmed ☐ Untrimmed

Spiced Alignment? ☒

Select Aligner

☐ Rbowtie ☒ Rhisat2

Max Mapping Positions Per Read

1

Read Position

☒ start ☐ end

Orientation of Alignments

☒ any ☐ same ☐ opposite

Read Mate to Count Alignments

☒ any ☐ first ☐ last

Count Spiced Alignments? ☒

Count Secondary Alignments? ☒

Mapping Quality of Alignments When Counting

0 255

0 26 52 78 104 130 156 182 208 234 255

Run Aligning

Click once done

Fig 4.16: All alignment settings

4. Click “Run Aligning” once done
(note: There’ll be 3 modal dialogs which appear: “Aligning”, “Generating Counts” and “Quality Reporting”)
5. In the “Results” tab panel, download the csv files for gene, exon, promoter and junction counts. You’ll also see a preview for these files.
(note: Junction counts are only generated when spliced alignment is enabled.)
(note: Column names are identifiers of the fastq files within the sample selected for aligning.)

Example of gene counts:

https://github.com/paigerollex/gene_cloud_omics/blob/main/output_data/male_a_gene_n_cbi.csv

Example of exon counts:

https://github.com/paigerollex/gene_cloud_omics/blob/main/output_data/male_a_exon_n_cbi.csv

Example of promoter counts:

https://github.com/paigerollex/gene_cloud_omics/blob/main/output_data/male_a_promoter_ncbi.csv

Example of junction counts:

https://github.com/paigerollex/gene_cloud_omics/blob/main/output_data/male_a_junction_ncbi.csv

The screenshot displays the user interface for uploading FASTQ files and viewing gene counts. On the left, a file explorer shows a folder named 'dm_male_a' containing three files: 'SRR999253.fastq.gz', 'SRR999254.fastq.gz', and 'SRR999255.fastq.gz'. The suffix '.fastq.gz' is highlighted in blue. On the right, a form for 'Sample Name' (male_a) and 'FASTQ Suffix' (.fastq.gz) is shown, along with a 'Select Single Read FASTQ File' section. Below the form, a 'Unzip FASTQ' button is visible. The main section, titled 'Gene Counts', includes a 'Download: Gene Counts' button and a 'Preview' table. The table has columns for 'width', 'SRR999253', 'SRR999254', and 'SRR999255'. The rows represent different genomic features: 128up, 14-3-3epsilon, 14-3-3zeta, 140up, 18SrRNA-Psi:CR41602, and 18SrRNA-Psi:CR45861.

	width	SRR999253	SRR999254	SRR999255
128up	1304	123	90	95
14-3-3epsilon	3366	6931	6213	6403
14-3-3zeta	3637	15361	14034	14520
140up	1200	100	74	79
18SrRNA-Psi:CR41602	1975	405	284	327
18SrRNA-Psi:CR45861	1934	3218	1841	2065

Fig 4.17: Gene counts

Exon Counts

Download: Exon Counts

Preview:

	width	SRR999253	SRR999254	SRR999255
1	201	0	1	0
10	588	2	0	2
100	166	58	58	73
1000	1713	683	637	708
10000	313	294	274	262
100000	142	3	6	1

Fig 4.18: Exon counts

Promoter Counts

Download: Promoter Counts

Preview:

	width	SRR999253	SRR999254	SRR999255
1;NM_001110622.3	2200	0	1	0
10;NM_001258476.2	2200	47	37	44
100;NM_001272138.1	2200	174	139	163
1000;NM_167020.2	2200	8	7	14
10000;NM_143232.3	2200	0	2	0
10001;NM_170281.2	2200	0	0	0

Fig 4.19: Promoter counts

Junction Counts

Download: Junction Counts

Preview:

	seqnames	start	end	width	strand	SRR999253	SRR999254
1	NC_004353.4	1004413	1004469	57	+	14	8
2	NC_004353.4	1004413	1004469	57	-	13	8
3	NC_004353.4	1004619	1005806	1188	+	19	21
4	NC_004353.4	1004619	1005806	1188	-	27	27
5	NC_004353.4	1006010	1006065	56	+	17	14
6	NC_004353.4	1006010	1006065	56	-	15	17

Fig 4.20: Junction counts

